

EXAMVISION™



Congratulations on your
custom-crafted loupe



Content

Introduction	3
Quality control, fine-adjustment and videos	
Galilean HD & Essential	4
Kepler Kompakt	6
Kepler Advanced	8
Warranty	10
ExamVision loupe accessories	11
Symbols and instructions for proper use	14

Introduction

You have chosen a high-quality precision tool tailored entirely to you. If your loupe needs adjusting, please remember it is a complex precision instrument, so follow the video guidance carefully – or ask your local ExamVision dealer for assistance.

ExamVision can adjust the working distance of your loupe, and upgrade prescription lenses as your vision alters over time. If you have any questions or concerns, please contact your ExamVision dealer.

Your loupe has been produced with care, to precisely your specifications, and we hope you are 100% satisfied.



Quality control, fine-adjustment and videos





Galilean HD & Essential

Before you start with your new loupe, please watch the ExamVision videos to ensure correct fitting and vision.

We have set your loupe to your desired working distance from the front of your eye to your subject. This distance is written on the Quality Card attached to your new loupe.

Watch our video "*Does your loupe fit correctly?*" for optimal clear vision, comfort and perfect fitting. This ensures the correct height and angle of the loupe, as illustrated below:



Videos

You can find ExamVision videos on our channel on YouTube and also on Vimeo. To watch now just scan the QR code here:





Quality control, fine-adjustment and videos



Kepler Kompakt

Before you start with your new loupe, please watch the ExamVision videos to ensure correct fitting and vision.

We have set your loupe to your desired working distance from the front of your eye to your subject. This distance is written on the Quality Card attached to your new loupe.

You must fine-focus your Kepler Kompakt one eye at a time, at your working distance. Fix your optimal focus setting using the enclosed Allen key. Watch our video *"Fitting your Kepler Kompakt loupe"* for optimal clear vision, comfort and perfect fitting. This ensures the correct height and angle of the loupe, as illustrated below:



Videos

You can find ExamVision videos on our channel on YouTube and also on Vimeo. To watch now just scan the QR code here:

 YouTube





Quality control, fine-adjustment and videos



Kepler Advanced

Before you start with your new loupe, please watch the ExamVision videos to ensure correct fitting and vision.

We have set your loupe to your desired working distance from the front of your eye to your subject. This distance is written on the Quality Card attached to your new loupe.

You must fine-focus your Kepler Advanced one eye at a time, at your working distance. Watch our video *"Fitting your Kepler Advanced loupe"* for optimal clear vision, comfort and perfect fitting. This ensures the correct height and angle of the loupe, as illustrated below:



Videos

You can find ExamVision videos on our channel on YouTube and also on Vimeo. To watch now just scan the QR code here:



Warranty

5 year manufacturer warranty

- All Galilean loupes
- Kepler Kompakt
- Kepler Advanced
- Essential

Normal wear and tear and damage due to inadequate use or maintenance are not covered by warranty. In the event of malfunction, please return the device to your ExamVision dealer for repair. Please also see examvision.com/warranty.

ExamVision loupe accessories



Nose pads

Clear silicone nose pads in various sizes. Included with all ExamVision loupes.



Neck cord

Delicate yet durable cord for extra security when wearing your loupes - with silicone attachments. Included with all ExamVision loupes.



Cable clip & adapters for visionshields

Attach the cable clips to the temples of your loupes and push your LED cable gently into the clip. Your visionshield can now be attached to your loupe. Included with all ExamVision loupes.



Cleaning solution (non-alcoholic)

Our custom cleaning solution will keep your loupe in perfect condition. Use a small amount on a microfibre lens cloth to safely remove dirt, dust and grease. Included with all ExamVision loupes.

ExamVision loupe accessories



Disinfectant wipes (non-alcoholic)

Our disinfectant wipes will keep your loupe in perfect condition. These disposable cleaning wipes are a quick and easy way to safely remove dirt, dust and grease.



Visionshield (size L, XL and MAX)

Protects the loupe and your face from splashes and dirt. Easily attached to ExamVision loupes using a Velcro adapter. Size L fits Galilean HD and Essential loupes – and size XL fits Kepler loupes. These versions have a cut out in the visor for a head lamp. The Visionshield MAX is a full-face visor is easily attached to and fits on any ExamVision loupe as part of your work routine. ExamVision visionshields are made from a non-splintering, flexible plastic film which can be easily cleaned, disinfected and reused.



Sideshields

Almost invisible, featherlight side-shields protect the eyes from splashes and debris. Disposable, in packs of three pairs.



Jewels

If total customisation is not enough, why not add a diamond? A 0.02 carat (1.7mm) diamond nestling in a sterling silver stud, handmade by a Danish jewellery designer to enhance your beautiful Icon frame. Also, rubies and sapphires are available. We perfect it. You personalise it.



Coloured temple tips

Anti-allergy temple tips for personalization and a comfortable fit. Made of artificial rubber. Available in a wide range of colours.

You can see the full list of available accessories on examvision.com/loupes-lights/accessories/

Symbols and instructions for proper use



Manufacturer

ExamVision ApS, Industrivej 11, 8305 Samsø, Denmark.
www.examination.com



Consult instructions for use or consult electronic instructions for use.
Read all enclosed guides and instructions carefully before use.

Find all user guides here:
examvision.com/support



Caution. Correct use and maintenance

Your ExamVision loupe is a precision instrument. As with any quality product, proper care is needed to keep your loupe in peak condition. With normal daily use, your loupe will need service checks from time to time. Inspect regularly, and if you notice wear or damage, if your loupe feels loose or does not fit as well as it did when new, discontinue use and contact your ExamVision dealer immediately for service and replacement parts. Your local eyecare professional may also be able to assist you with small adjustments or tightening of screws. Do NOT modify or alter your loupe.

Allergies

ExamVision uses only high-quality, allergy-tested eyeglass materials. Your loupe must be fitted so only the nose pad and temple tips touch your skin. In the highly unlikely event of any irritation, contact your ExamVision dealer.

Care and cleaning

Do NOT use alcohol for cleaning. Use only ExamVision recommended alcohol-free cleaning and disinfectant products and a soft lint-free cloth. Do NOT autoclave or use ultra-sonic cleaning methods. Do NOT use UV immersive disinfective methods.



Keep away from sunlight

Do not leave loupe glasses directly in the sun. Magnification can create a "hot spot" which can cause a fire.



Keep dry.



Fragile, handle with care.



Do not use if package is damaged.



 facebook.com/ExamVisionDK

 instagram.com/examvision/

www.examvision.com