

# Focus Bright & Xtend

## Technical Specifications



Product name:	Focus Bright & Xtend
Weight	Lamp: 17 g. Battery: 156 g.
Light colour	Bright colour, neutral white ( $\leq 6.500$ kelvin)
Brightness settings	Low   Medium   High
Light intensity at 35 cm	18,600 lux
Light spot at 35 cm	70 mm.
Heat development	40 to 60°C, Depending on setting and time in use
Operating time	6-24 hours depending on chosen light intensity
Charging time	Approx. 2.5 hours
Battery type	Focus Xtend Control Unit
Battery pack	Lithium-ion, 2 cell 7,2V 1.2 AH
Power supply	12-13 V DC, 1A Medico
Temperature range	Operating: 10 – 35°C Storage: -20 – 60°C Humidity: 10 – 95 %
Warranty	3 years on Control Unit Xtend, Battery Xtend, lamps (Focus True, Focus Bright, & Focus Nano) and accessories. 1 year on detachable cable and backplate