

Focus Bright & Xtend

Technical Specifications



Product name:	Focus Bright & Xtend
Weight	Lamp: 17 g. Battery: 156 g.
Brightness settings	Low Medium High
Operating time	6-24 hours depending on chosen light intensity
Charging time	Approx. 2.5 hours
Light intensity at 35 cm	18,600 lux
Light spot at 35 cm	70 mm.
Light colour	Bright colour, neutral white (≤ 6.500 kelvin)
Heat development	40 to 60°C, Depending on setting and time in use
Battery type	Focus Xtend Control Unit
Battery pack	Lithium-ion, 2 cell 7,2V 1.2 AH
Power supply	12-13 V DC, 1A Medico
Temperature range	Operating: 10 – 35°C Storage: -20 – 60°C Humidity: 10 – 95 %
Warranty	3 years on Control Unit Xtend, Battery Xtend, lamps (Focus True, Focus Bright, & Focus Nano) and accessories. 1 year on detachable cable and backplate