

Please note that the actual perceived magnification can vary, depending on working distance and prescription.

These examples are to be used as a guide and represent angular magnification, based on a working distance of 45 cm.



No magnification

## Galilean HD



<b>Magnification</b>	<b>2.3 x</b>	<b>2.8 x</b>	<b>3.3 x</b>	<b>3.8 x</b>
Field of View	12 cm	10.5 cm	8 cm	6.5 cm
Depth of Field up to	28 cm	20 cm	15 cm	10 cm

## Essential



<b>Magnification</b>	<b>2.5 x</b>	<b>3.0 x</b>
Field of View	11 cm	8.3 cm
Depth of Field up to	25 cm	15 cm