

Please note that the actual perceived magnification can vary, depending on working distance and prescription.

These examples are to be used as a guide and represent angular magnification, based on a working distance of 45 cm.



No magnification

Galilean HD



Magnification	2.3x	2.8x	3.3x	3.8x
Field of View	12 cm	10.5 cm	8 cm	6.5 cm
Depth of Field up to	28 cm	20 cm	15 cm	10 cm

Galilean Essential



Magnification	2.5x	3.0x
Field of View	11 cm	8.3 cm
Depth of Field up to	25 cm	15 cm