



Kepler Reflekt

Technical Specifications

Loupe lens system	Kepler Reflekt
Lenses	Plano or Prescription
BlueLight Protection (BLP)	Yes (optional add-on)
Prisma Lens System	Yes (optional add-on)
Fine-focus available	Yes
Magnifications	4.0x 6.0x
Working distance	4.0x: From 35 to 57 cm 6.0x: From 47 to 69 cm
Field of View (at 45 cm)	4.0x: 7.2 cm 6.0x: 5.6 cm
Depth of field (up to)	4.0x: 13 cm 6.0x: 8 cm
Mount system	TTL
Weight of two oculars with correction lenses	4.0x: 36.9 g 6.0x: 39.8 g
Loupe colour	Black housing with titanium front
Loupe material	Magnesium and titanium
Frames available	All