

Essential Alpha LED & Essential Alpha Control Unit

Technical Specifications



Product name:	Essential Alpha
Material	Anodised Aluminium
Weight	11 g
Brightness settings	Stepless
Operating time	Up to 12 hours
Charging time	Approx. 8 hours
Light intensity	At 25 cm: 37,000 lux At 35 cm: 14,000 lux
Light spot at 35 cm	85 mm
Light colour	Natural white (6,300 K)
Charger output	DC 5V.1A.5W
Cable type	Black super flex
Warranty	Battery and cable 1 year Lamp 2 years