



Kepler Advanced Technical Specifications

Variant	Kepler Advanced v2 – 4.0x + 5.0x + 6.0x + 7.0x
Lenses	Plano or Prescription
BlueLight Protection (BLP)	Yes (optional add-on)
Prisma Lens System	Yes (optional add-on)
Fine-focus available	Yes
Magnification	Turn’N’Click option of: 4.0x + 5.0x + 6.0x + 7.0x
Working distance	From 31 cm to 69 cm
Field of view (at 45 cm)	4.0x: 6.5 cm 5.0x: 6.0 cm 6.0x: 5.5 cm 7.0x: 5.0 cm
Depth of field (up to)	4.0x: 8 cm 5.0x: 7 cm 6.0x: 6 cm 7.0x: 5 cm
Mount System	TTL (Through the lens)
Weight of two oculars with correction lenses	49.5 g
Loupe colour	Anthracite grey/Metallic grey
Loupe material	Carbon fibre reinforced composite and Titanium
Frames available	All