

EXAM **V**ISION

New

Kepler Reflekt

Your vision is
our mission




Unobstructed view

No need to remove the loupes to perform other tasks.

Superior ergonomics & peripheral vision

Kepler Reflekt is the first reflective ergo loupe by ExamVision. The slim black loupe housing is made from durable yet lightweight magnesium and titanium.

Unique to its streamlined ergonomic design is the unmatched peripheral vision offering full, clear visibility of your surroundings and allowing you to continue your daily work without the need to remove the loupes.



Peripheral vision obtained with Kepler Reflekt compared to conventional periscope loupes.



We flipped the prism

This allows you to see above the ocular without obstruction.

Lens customisation

As with all ExamVision loupes we offer the finest lenses compatible with all needs.

Magnesium and titanium oculars

Robust yet lightweight in an ultra-narrow design making it even lighter.

Adjustable fine-focus

A standard on ExamVision Kepler loupes which gives you that extra level of customisation by allowing each ocular to be set perfectly to the individual sight of each eye.

4.0x and 6.0x magnification

Suitable for prolonged use and complex procedures.

Natural and comfortable ergonomics





The Kepler Reflekt allows you to sit with a lower neck angle than conventional direct-vision loupes, putting less weight on your neck throughout the day.





Your vision
is our mission

EXAMVISION

-  facebook.com/examvisiondk
-  instagram.com/examvision/
-  linkedin.com/company/examvision-aps
-  youtube.com/@examvisionchannel

www.examvision.com

Scan the QR code for
more information

