

Kepler Kompakt and Kepler Advanced are direct view (prismatic) Kepler systems. These examples are to be used as a guide and represent a linear magnification, based on a working distance of 45 cm.



No magnification

Kepler Kompakt



Magnification

3.4x

5.0x

7.6x

Field of View

9.5 cm

6.2 cm

4.3 cm

Depth of Field up to

11.5 cm

8 cm

4.3 cm

Kepler Advanced



Magnification

4.0x

5.0x

6.0x

7.0x

Field of View

6.5 cm

6 cm

5.5 cm

5 cm

Depth of Field up to

8 cm

7 cm

6 cm

5 cm