



PowerGo

Technical Specifications

Control Unit	PowerGo
Material	Aluminium
Weight	155 g
Compatible LED Lights	Total Pure & Total Intense LED Lights
Connection type	USB-C
Light intensity settings	High, Medium, and low
Operating time	Minimum 6 hours at high intensity
Charging time	Approx. 2-2.5 hours
Battery Type	21,700 Lithium-ion rechargeable battery
Output rating	3 VDC 0.5-2.0 A
Input rating	5 VDC 2.0-2.5 A
Warranty	PowerGo Unit: 3 years PowerGo Battery: 1 year