



Galilean Essential

Technical Specifications

Loupe lens system	Galilean Essential
Lenses	Plano or Prescription
BlueLight Protection (BLP)	Yes (optional add-on)
Prisma Lens System™	Yes (optional add-on)
Magnification	2.5x 3.0x
Working distance	From 30 to 55 cm
Field of view (at 45 cm)	2.5x: 9.5 cm 3.0x: 7.5 cm
Depth of field (at age 40 years)	2.5x: 25 cm 3.0x: 15 cm
Mount system	TTL (Through the Lens)
Weight (in a Panto Small frame)	2.5x: 38.6 g 3.0x: 40.5 g
Loupe colour	Metallic Grey
Loupe material	Ingeo™ Biopolymer
Frames available	All