



Galilean Essential

Technical Specifications

Loupe lens system	Galilean Essential
Lenses	Plano or Prescription
BlueLight Protection (BLP)	Yes (optional add-on)
Prisma Lens System	Yes (optional add-on)
Fine-focus available	No
Magnification	2.5x 3.0x
Working distance	From 30 to 55 cm
Field of view (at 45 cm)	2.5x: 11 cm 3.0x: 8.3 cm
Depth of field (up to)	2.5x: 25 cm 3.0x: 15 cm
Mount system	TTL (Through the Lens)
Weight of two oculars with correction lenses	2.5x: 19.4 g 3.0x: 21.4 g
Loupe colour	Metallic grey
Loupe material	Ingeo™ Biopolymer
Frames available	All