



Galilean Essential Technical Specifications

Variant	Galilean Essential – 2.5x & 3.0x	
Lenses	Plano or Prescription	
BlueLight Protection (BLP)	Yes (optional add-on)	
Prisma Lens System	Yes (optional add-on)	
Fine-focus available	No	
Magnification	2.5x	3.0x
Working distance	From 30 cm to 55 cm	
Field of view (at 45 cm)	2.5x: 11 cm	3.0x: 8.3 cm
Depth of field (up to)	2.5x: 25 cm	3.0x: 15 cm
Mount System	TTL (Through the lens)	
Weight of two oculars with correction lenses	2.5x: 19.4 g	3.0x: 21.4 g
Loupe colour	Metallic grey	
Loupe material	Ingeo Biopolymer	
Frames available	All	