



Essential

Technical Specs

Loupe lens system	Galilean
Lenses	Plano or Prescription
BlueLight Protection (BLP)	Yes (optional add-on)
Prisma Lens System™	Yes (optional add-on)
Magnification	2.5x 3.0x
Working distance	2.5x From 33 to 60 cm 3.0x From 35 to 50 cm
Field of view (at 40 cm)	2.5x: 8.7 cm 3.0x: 5.5 cm
Depth of Field (at age 40 years)	2.5x: 18 cm 3.0x: 14 cm
Mount system	TTL
Weight (in a Titanium Cosmo frame, Small)	2.5x: 52 g 3.0x: 56 g
Weight (in a Monel Essential frame, Small)	2.5x: 62 g 3.0x: 65 g
Loupe colour	Grey
Loupe material	Aluminum housing
Frames available	All