



# Kepler Kompakt

## Technical Specifications

Loupe lens system	Kepler Kompakt
Lenses	Plano or Prescription
BlueLight Protection (BLP)	Yes (optional add-on)
Prisma Lens System	Yes (optional add-on)
Fine-focus available	Yes
Power-Up System	Yes, front ocular is changeable
Magnification	3.5x   4.6x   5.7x
Working distance	From 30 to 60 cm
Field of view (at 45 cm)	3.5x: 8.5 cm   4.6x: 7.5 cm   5.7x: 6.5 cm
Depth of field (up to 50 years)	3.5x: 10 cm   4.6x: 8 cm   5.7x: 6 cm
Mount system	TTL (Through the Lens)
Weight of two oculars	3.5x: 57.8 g   4.6x: 57 g   5.7x: 58.6 g
Loupe colour	Titanium grey housing, TTL black composite
Loupe material	Titanium and composite housing
Frames available	All (except Carbon)