



Kepler Kompakt

Technical Specifications

Loupe lens system	Kepler Kompakt
Lenses	Plano or Prescription
BlueLight Protection (BLP)	Yes (optional add-on)
Prisma Lens System	Yes (optional add-on)
Fine-focus available	Yes
Power-Up System	Yes, front ocular is changeable
Magnification	3.5x 4.6x 5.7x
Working distance	From 30 to 60 cm
Field of view (at 45 cm)	3.5x: 8 cm 4.6x: 7.5 cm 5.7x: 5.5 cm
Depth of field (up to)	3.5x: 9 cm 4.6x: 6.5 cm 5.7x: 6 cm
Mount system	TTL (Through the Lens)
Weight of two oculars with correction lenses	3.5x: 57.8 g 4.6x: 57 g 5.7x: 58.6 g
Loupe colour	Titanium and black
Loupe material	Titanium and black anodised aluminium
Frames available	All (except Carbon)