

EXAM **V**ISION




Extraordinary performance

Finest, custom-crafted Danish loupe
and LED systems



HANDMADE
in Denmark



Angulation tailored to your primary needs, whether that involves mainly small animals, large animals or surgery



Danish design. Handcrafted for veterinary professionals.

At ExamVision, Danish craftsmanship and optical skill meet veterinary expertise. Every loupe is handmade in Denmark using premium materials, such as titanium, carbon fibre and stainless steel, and individually customised to your prescription, working distance and ergonomics.

ExamVision means:

- Visual clarity & ergonomic improvement
- Handcrafted Danish quality with premium materials
- Individually customised for vision, posture and prescription
- Rebuildable to follow prescription changes

Veterinary benefits:

- More detail at a shorter, safer working distance for accurate diagnostics and treatment
- Galilean HD or Galilean Essential loupes, both compact, robust, lightweight and suitable for all-day wear
- Alpha LED light with bright, natural illumination and adjustable intensity for every procedure



Try our 3D Loupe Creator and bring your dream loupe to life in 4 simple steps.



Galilean Essential

- The ideal entry-level loupe offering all-day comfort and precision
- High-quality magnification at an attractive price
- Lightweight thanks to shorter ocular housings made from Ingeo™ biopolymer
- Available magnifications: 2.5x, 3.0x

Individual magnifications

2.5 3.0



Galilean HD

- Perfectly balanced comfort for optimal use and posture
- Maximum depth of field and exceptionally wide field of view
- Ideal for all-day wear
- Frames and oculars crafted in titanium
- Adaptable to any prescription
- Prisma Lens System available
- Available magnifications: 2.3x, 2.8x, 3.3x, 3.8x

Individual magnifications

2.3 2.8 3.3 3.8



Essential Alpha LED

- Natural white light $\leq 6,300$ Kelvin
- Light intensity 14,000 Lux at 35 cm
- Light spot diameter 85 mm at 35 cm
- Continuously adjustable brightness
- Lamp weight only 11 g
- Up to 12 hours of operating time



“
*ExamVision loupes
are so tailor made
and comfortable,
that I use them for
every single consult*

Christian Fabricius
CMV, Veterinarian, Denmark

”

Learn more about
our magnifications



Contact your local
ExamVision dealer



EXAMVISION

-  facebook.com/examvisiondk
-  instagram.com/examvision
-  linkedin.com/company/examvision-aps
-  youtube.com/@examvisionchannel

www.examvision.com



25.11.03:EN