

Essential Nano & Focus Xtend

Technical Specifications



Product name:	Essential Nano
Weight	Lamp: 9 g. Battery: 156 g.
Brightness settings	Variable
Operating time	Up to 12 hours depending on chosen light intensity
Charging time	Approx. 2.5 hours
Light intensity at 35 cm	27,400 lux
Light spot at 35 cm	80 mm.
Light colour	Natural white (≤ 6.500 kelvin)
Heat development	40 to 60°C, Depending on setting and time in use
Battery type	Focus Xtend Control Unit
Battery pack	Lithium-ion, 2 cell 7,2V 1.2 AH
Power supply	12-13 V DC, 1A Medico
Temperature range	Operating: 10 – 35°C Storage: -20 – 60°C Humidity: 10 – 95 %
Warranty	3 years on Control Unit Xtend, Battery Xtend, lamps (Focus True, Focus Bright, & Focus Nano) and accessories. 1 year on detachable cable and backplate