



Extraordinary performance

Finest, custom-crafted Danish loupe and LED systems





Custom-crafted medical loupes for your needs

ExamVision offers versatile, all-day-wear tools to physicians, surgeons, and other medical professionals and students who require the ultimate optical loupe. Remarkably lightweight yet made from the most robust, high-end materials available.

Medical examinations and procedures require the use of superior optics and tools to facilitate freedom of movement, good posture and correct ergonomics throughout your working day.

Choosing the right loupe and magnification will optimise your performance in surgical and other medical procedures. Each solution is based on individual eyesight, need for magnification as well as physical working environment such as distance, tools and facilities available. So, your loupe is truly tailored and custom-crafted to you.

Your benefits

All of our products are the result of a perfect synergy between technology and handcraft. We produce and assemble each loupe based on individual measurements. They then undergo a rigorous quality control process before leaving the factory.

- Excellent vision and precision
 Our optics help you achieve the crystal-clear image you need for high-quality, precision work.
- Work satisfaction
 Achieved through better treatment results, faster work and happier, well-treated patients.
- Peace of mind
 Five-year warranty on all of our loupes, with unrivalled dealer service support.

- Longer working life
 An improved ergonomic position via customised loupes is proven to avoid musculoskeletal problems.
- Style and functionality
 Six-time winner of the international product design award Red Dot Design Award.
- Customised solution
 Our specialists in optics and ergonomics ensure a solution tailored uniquely to you.

ExamVision loupe systems



Kepler Advanced

- The only loupe with 4 magnifications
- · Adjustable fine-focus
- · Switch mag. easily via the Turn'N'Click system
- · Perfect magnification level for all your demanding procedures
- Unparalleled sharpness from low to high magnification
- Made of titanium and carbon fibre reinforced composite

4 magnifications in 1 loupe









Kepler Reflekt

- · Magnesium and titanium materials ensure a robust and ultra-lightweight loupe
- · Improved ergonomics thanks to steep angulation
- · Wide field of view, crystal-clear from edge to edge
- 4.0x magnification for all-day wear and detailed procedures
- · Adjustable fine-focus for highest level of personalisation

Individual magnification







- Perfectly balanced comfort ensuring optimum use and posture
- · Maximum depth of field
- · Impressively wide field of view allowing plenty of movement
- · Compact and extremely lightweight - ideal for all day wear
- · Titanium oculars available in three colours matching your frame colour

Individual magnifications











Kepler Kompakt

- · Lightweight, ultra-compact design
- · Wide field of view, crystal clear from edge to edge
- PowerUp system to upgrade or downgrade magnification level
- Lockable fine-focus ensures precision
- Freedom of movement maintained by excellent depth of field

Individual magnifications









ExamVision Total LEDs

Having the right light at the right spot is vital both in terms of good vision and protecting your eyes from glare and harmful light. Loupes and light are interdependent and choosing the right headlight solution is almost as vital as choosing the optimum eye loupe solution. Your ExamVision dealer will guide you through personal consultation to ensure that you get the maximum benefit from your system.



Total Intense

- Bright, natural white 6,500 Kelvin
- CRI above 70
- Light intensity >25,700 lux @ 35 cm
- Light spot, 70 mm @ 35 cm
- Three brightness settings
- Lamp weight 18 g

Total Pure

- Bright, natural white 5,000 Kelvin
- CRI above 90
- Light intensity >11,700 lux @ 35 cm
- Light spot, 70 mm @ 35 cm
- Three brightness settings
- Lamp weight 18 g



The Total LED comes with the option of a coaxial attachment. This will ensure the light follows your line of sight if you need a steep loupe angle.

Ergonomics

Good ergonomics are essential for a long, enjoyable, healthy working life.

Ergonomics and safety

By using magnification, you will not only achieve an optimal, ergonomic working position but also increase the safety distance between you and your patient significantly. Magnifying loupes will allow you to perform professionally without compromising your own body or protection.

Get qualified ergonomics advice

With your ExamVision loupes, you also receive a free session with your specialist dealer. During this session, your dealer offers you ergonomical advice to help you work in a way that promotes physical well-being and good posture.

Prisma Lens System

Prisma Lens System is optional in your oculars. It allows you to sit up straighter and achieve a better working posture by making the focus of your vision appear higher than it really is.



Try our 3D loupe creator



Learn more about our magnifications



Contact your local ExamVision dealer





- f facebook.com/examvisiondk
- (instagram.com/examvision
- in linkedin.com/company/examvision-aps
- youtube.com/@examvisionchannel