

EXAM **V**ISION



Extraordinary performance

Finest, custom-crafted  
Danish loupe and LED systems

 **HANDMADE**  
in Denmark



### Our essence

At ExamVision we combine the rich heritage of Danish design with a commitment both to aesthetics and functionality. We come from Samsø, a tiny Danish island which is known worldwide for its renewable energy. The peace and natural beauty of Samsø affords us the perfect environment in which to think big, work hard and pay meticulous attention to every single detail.



## Custom-crafted medical loupes for your needs

ExamVision offers versatile, all-day-wear tools to physicians, surgeons, and other medical professionals and students who require the ultimate optical loupe. Remarkably lightweight yet made from the most robust, high-end materials available.

Medical examinations and procedures require the use of superior optics and tools to facilitate freedom of movement, good posture and correct ergonomics throughout your working day.

Choosing the right loupe and magnification will optimise your performance in surgical and other medical procedures. Each solution is based on individual eyesight, need for magnification as well as physical working environment such as distance, tools and facilities available. So, your loupe is truly tailored and custom-crafted to you.

## Your benefits

All of our products are the result of a perfect synergy between technology and handcraft. We produce and assemble each loupe based on individual measurements. They then undergo a rigorous quality control process before leaving the factory.

- **Excellent vision and precision**  
Our optics help you achieve the crystal-clear image you need for high-quality, precision work.
- **Longer working life**  
An improved ergonomic position via customised loupes is proven to avoid musculoskeletal problems.
- **Work satisfaction**  
Achieved through better treatment results, faster work and happier, well-treated patients.
- **Style and functionality**  
Six-time winner of the international product design award – Red Dot Design Award.
- **Peace of mind**  
Five-year warranty on all of our loupes, with unrivalled dealer service support.
- **Customised solution**  
Our specialists in optics and ergonomics ensure a solution tailored uniquely to you.

# ExamVision loupe systems



NEW

## Kepler Advanced

- The only loupe with 4 magnifications
- Adjustable fine-focus
- Improved working distance for better ergonomics
- Versatility across all procedures
- Unparalleled sharpness from low to high magnification
- Handcrafted from premium titanium and carbon fibre

4 magnifications in 1 loupe

4.0 5.0 6.0 7.0



NEW

## Kepler Kompakt

- Lightweight, ultra-compact design
- Offers the widest magnification range in the ExamVision lineup – up to 7.6x
- Wide field of view, crystal-clear from edge to edge
- Adjustable fine-focus for perfect clarity
- Excellent depth of field allows natural freedom of movement

Individual magnifications

3.4 5.0 7.6

NEW

## Kepler Reflekt

- Magnesium and titanium materials ensure a robust and ultra-lightweight loupe
- Improved ergonomics thanks to steep angulation
- Wide field of view, crystal-clear from edge-to-edge
- 4.0x and 6.0x magnification for all-day wear and detailed procedures
- Adjustable fine-focus for perfect clarity

Individual magnifications

4.0 6.0



reddot winner 2024

## Galilean HD

- Perfectly balanced comfort ensuring optimum use and posture
- Maximum depth of field
- Impressively wide field of view allowing plenty of movement
- Compact and extremely lightweight – ideal for all day-wear
- Titanium oculars available in three colours matching your frame colour

Individual magnifications

2.3 2.8 3.3 3.8



reddot winner 2010



## ExamVision Total LEDs

Having the right light at the right spot is vital both in terms of good vision and protecting your eyes from glare and harmful light. Loupes and light are interdependent and choosing the right headlight solution is almost as vital as choosing the optimum eye loupe solution. Your ExamVision dealer will guide you through personal consultation to ensure that you get the maximum benefit from your system.



### Total Intense

- Bright, natural white 6,500 Kelvin
- CRI above 70
- Light intensity >25,700 lux @ 35 cm
- Light spot, 70 mm @ 35 cm
- Three brightness settings
- Lamp weight 18 g

### Total Pure

- Bright, natural white 5,000 Kelvin
- CRI above 90
- Light intensity >11,700 lux @ 35 cm
- Light spot, 70 mm @ 35 cm
- Three brightness settings
- Lamp weight 18 g



The Total LED comes with the option of a coaxial attachment. This will ensure the light follows your line of sight if you need a steep loupe angle.

## Ergonomics

Good ergonomics are essential for a long, enjoyable, healthy working life.

### Ergonomics and safety

By using magnification, you will not only achieve an optimal, ergonomic working position but also increase the safety distance between you and your patient significantly. Magnifying loupes will allow you to perform professionally without compromising your own body or protection.

### Get qualified ergonomics advice

With your ExamVision loupes, you also receive a free session with your specialist dealer. During this session, your dealer offers you ergonomical advice to help you work in a way that promotes physical well-being and good posture.

### Prisma Lens System

Prisma Lens System is optional in your oculars. It allows you to sit up straighter and achieve a better working posture by making the focus of your vision appear higher than it really is.




Try our 3D  
loupe creator



Learn more about  
our magnifications



Contact your local  
ExamVision dealer



-  [facebook.com/examvisiondk](https://facebook.com/examvisiondk)
-  [instagram.com/examvision](https://instagram.com/examvision)
-  [linkedin.com/company/examvision-as](https://linkedin.com/company/examvision-as)
-  [youtube.com/@examvisionchannel](https://youtube.com/@examvisionchannel)

[www.examvision.com](https://www.examvision.com)

