



Galilean HD

Technical Specifications

Loupe lens system	Galilean HD
Lenses	Plano or Prescription
BlueLight Protection (BLP)	Yes (optional add-on)
Prisma Lens System	Yes (optional add-on)
Fine-focus available	No
Magnification	2.3x 2.8x 3.3x 3.8x
Working distance	From 30 to 60 cm
Field of view (at 45 cm)	2.3x: 12 cm 2.8x: 10.5 cm 3.3x: 8 cm 3.8x: 6.5 cm
Depth of field (up to)	2.3x: 28 cm 2.8x: 20 cm 3.3x: 15 cm 3.8x: 10 cm
Mount system	TTL (Through the Lens)
Weight of two oculars with correction lenses	2.3x: 31.8 g 2.8x: 26.2 g 3.3x: 34.4 g 3.8x: 34.4 g
Loupe colour	Titanium grey
Loupe material	Titanium housing
Frames available	All