



Galilean HD

Technical Specifications

| | |
|--|---|
| Variant | Galilean HD |
| Lenses | Plano or Prescription |
| BlueLight Protection (BLP) | Yes (optional add-on) |
| Prisma Lens System | Yes (optional add-on) |
| Fine-focus available | No |
| Magnification | 2.3x 2.8x 3.3x 3.8x |
| Working distance | From 30 cm to 60 cm |
| Field of view (at 45 cm) | 2.3x: 12 cm 2.8x: 10.5 cm 3.3x: 8 cm 3.8x: 6.5 cm |
| Depth of field (up to) | 2.3x: 28 cm 2.8x: 20 cm 3.3x: 15 cm 3.8x: 10 cm |
| Mount System | TTL (Through the lens) |
| Weight of two oculars with correction lenses | 2.3x: 31.8 g 2.8x: 26.2 g 3.3x: 34.4 g 3.8x: 34.4 g |
| Loupe colour | Titanium Grey, Rose Gold & Blue Slate |
| Loupe material | Titanium |
| Frames available | All |