



# Galilean HD

## Technical Specs

Loupe lens system	Galilean HD
Lenses	Plano or Prescription
BlueLight Protection (BLP)	Yes (optional add-on)
Prisma Lens System™	Yes (optional add-on)
Fine adjustable	No
Lockable Fine-focus available	No
Magnification	2.3x   2.8x   3.3x   3.8x
Working distance	From 30 to 60 cm
Field of view (at 40 cm)	2.3x: 10.8 cm   2.8x: 10.2 cm   3.3x: 8.13 cm   3.8x: 7.11 cm
Depth of Field (at age 50 years)	2.3x: 20 cm   2.8x: 15 cm   3.3x: 10 cm   3.8x: 5 cm
Mount system	TTL
Weight (in an Icon Small frame)	2.3x & 2.8x, 2 layers lenses = 55 g 3.3x & 3.8x, 3 layers lenses = 60 g
Loupe colour	Titanium grey
Loupe material	Titanium housing
Frames available	All